



## Genotyping PCR for SIK3 KO mice

### 1. Method used in our animal resource bank

This section describes our method we are routinely using for genotyping SIK3 KO mice.

#### 1.1) Primer sequences:

- Primer1:a
  - Sequence: 5'- GCT ACC AAC TTG GTT ACA GTT GCT -3' (24-mer)
- Primer2: b2
  - Sequence: 5'- CGA GGG TCA GCA GCA ACT TCT-3' (21-mer)
- Primer3: c2
  - Sequence: 5'- GAG ACT AGT GAG ACG TGC TAC T-3' (22-mer)

#### 1.2) Reaction mixture:

	Tube (μL)
Water	8
Primer1 (a,10 μM)	0.5
Primer2 (b2,10 μM)	0.25
Primer3 (c2,10 μM)	0.25
Taq polymerase	10
DNA extracted from tail (diluted 500 times)	1
total	20

Taq polymerase: HotStarTaq Master Mix Kit (Qiagen). The enzyme is a chemically modified Taq polymerase for hot start PCR and needs 15-min incubation at 95 °C for activation. Master Mix contains enzyme, dNTP, Mg, etc at 2 x concentration. Please see Qiagen's website for details

(<http://www1.qiagen.com/Products/Pcr/HotStarTaqSystem/HotStarTaqMasterMix.aspx>).

#### 1.3) Thermal conditions:

94 °C	2 min	Enzyme activation and first denature
94 °C	30 sec	40 cycles
60 °C	30 sec	
72 °C	40 sec	
72 °C	2.5 min	once
4 °C	∞	

Thermal cycler: Veriti (Life technologies) with 0.2mL tubes.

#### 1.4) Product size:

- 250 bp for wild-type alleles
- 450 bp for KO alleles

Nibio\_ID=157;02/12/2013