



Genotyping of TNF KO mice

1. Method used in our animal resource bank

This section describes our method we are routinely using for genotyping TNF KO mice.

1.1) Primer sequences:

- Primer1: TNF Sense1
 - Sequence: 5'- AGA TGG AGA AGG GCA GTT AG -3' (20-mer)
- Primer2: TNF AntiS1
 - Sequence: 5'- ATA CCA GGG TTT GAG CTC AG -3' (20-mer)
- Primer3: TNF AntiS2
 - Sequence: 5'- TAC TTT GTT AAG AAG GGT GAG A -3' (22-mer)

1.2) Reaction mixture:

| | For wild allele | For KO allele |
|------------------------------------|-----------------|---------------|
| | Tube 1 (μL) | Tube 2 (μL) |
| Water | 8 | 8 |
| Primer1 (TNF Sense1, 10 μM) | 0.5 | 0.5 |
| Primer2 (TNF AntiS1, 10 μM) | 0.5 | - |
| Primer3 (TNF AntiS2, 10 μM) | - | 0.5 |
| Taq polymerase (U/μL) | 10 | 10 |
| DNA extracted from tail (purified) | 1 | 1 |
| total | 20 | 20 |

Taq polymerase: HotStarTaq Master Mix Kit (Qiagen). The enzyme is a chemically modified Taq polymerase for hot start PCR and needs 15-min incubation at 95 °C for activation. Master Mix contains enzyme, dNTP, Mg, etc at 2 x concentration. Please see Qiagen's website for details (<http://www1.qiagen.com/Products/Pcr/HotStarTaqSystem/HotStarTaqMasterMix.aspx>).

1.3) Thermal cycles:

Tube 1 (WT : primer1 and 2)

| | | |
|-------|--------|--------------------------------------|
| 95 °C | 15 min | Enzyme activation and first denature |
| 94 °C | 30 sec | |
| 60 °C | 1 min | 35 cycles |
| 72 °C | 1 min | |
| 72 °C | 10 min | |
| 4 °C | ∞ | once |

Tube 2 (KO : primer1 and 3)

| | | |
|-------|--------|--------------------------------------|
| 95 °C | 15 min | Enzyme activation and first denature |
| 94 °C | 30 sec | 35 cycles |
| 56 °C | 1 min | |
| 72 °C | 1 min | |
| 72 °C | 10 min | |
| 4 °C | ∞ | |

Thermal cycler: Veriti with 0.2mL tubes.

1.4) Product size:

Primers 1 and 2: 800 bp for wild-type alleles

Primers 1 and 3: 500 bp for KO alleles




Fig.1. Electropherogram of PCR products from wild and heterozygous mice.

Nibio_ID=156;11/07/2012